

## INSTRUCTIONS FOR USE

1. Inspect face mask and adapter for integrity of all parts before wearing the face mask. Examine for any gross contamination or the presence of foreign bodies. Breaks/cracks/tears can lead to ineffective protection against environmental contamination.



2. Before trying the face mask on for the first time, loosen all four black straps where they attach to the face mask (two connections at the top and two connections at the bottom).



3. Unclip the head strap clip by pushing the center button and gently pushing the clip apart.



4. Gently place the face mask against your face and align the face skirt of the face mask to your face.



5. Once the face mask feels comfortable on your face grab the two head strap clips and clip them together. You can have someone help you with this. We encourage practice with this technique until the user is comfortable with donning the face mask.



6. Ensure the clip is comfortable and seated on the back of your head.



7. Adjust and tighten the four head strap face mask connections so all four straps are pulling evenly.

